

## MCC – Milestone Cryoembedding Compound

### 1. Designated uses

Milestone Cryoembedding Compound is an excellent choice for supporting tissue during cryotomy.

Embedding medium for frozen tissue specimen to ensure Optimal cutting temperature.

Use Milestone Cryoembedding Compound to bind tissue to specimen block and to surround and encase specimen.

This is a well know water soluble glycol's and resins compound that provides an excellent specimen matrix for cryostat sectioning.

For in vitro diagnostic.

### 2. Waste disposal

For further information, please refer to the MSDS (Material Data Safety Sheet).

### 3. Symbol used



An instruction accompanied by this symbol provides a quite important and require more attention.



CE logo: this instrument complies to European Community directives.



IVD medical device.



Reference (code).



Lot number.



Manufacturer.



Expiration date.

